

MULTIPURPOSE OPERATION TABLE MOUNTED RETRACTOR SYSTEM



Serrated Slide-on Joint




TABLE POST



Clip-On Quick Angle 8"



Micro-Adjustable Quick Angle 10"



Malleable 62mm x 259mm (2 ½" x 10") S



Balfour [with lips] 30mm x 30mm (1" x 1")



Balfour [small lips] 71mm x 75mm (2 ¾" x 3")



Kelly 77mm x 88mm (3" x 3 ½")



Richardson 50mm x 175mm (2" x 7")



Kelly 62mm x 74mm (2 ½" x 3")



Balfour [with lips] 82mm x 72mm (3 ¼" x 3")



Non-Slip Balfour Right 80mm x 100mm (3" x 3 ¾")



Balfour [small lips] 76mm x 72mm (3" x 3")



Balfour [with lips] 81mm x 124mm (3 ¼" x 5")



Malleable 15mm x 210mm (½" x 8")



Malleable Finger 154mm (6")



Harrington 66mm x 155mm (2 ½" x 6")



Wide Balfour [no lips] 84mm x 58mm (3¼" x 2¼")



Wide Balfour [no lips] 110mm x 69mm (4" x 2 ⅜")



TABLE POST WITH CAM JOINTS

BILATERAL CROSS BAR STANDARD

BILATERAL CROSS BAR HINGED (ADJUSTABLE)

ANGLED SIDE ARM

SNAP LATCH CLIP-ON

SNAP LATCH WITH GEAR RACHET TYPE

RICHARDSONBLADE

BALFOUR BLADE

KELLY BLADE

HARRINGTON BLADE

ST.MARK'S PELVIC BLADE

MALLEABLE BLADE

FENESTRATED BLADE

DEAVER BLADE

(CUSTOMIZED BLADES) ON ORDER